



TIME CLOCK TUTORIAL

HOW TO USE THE TIME CLOCK

Use this feature if you want your loggers switched off when equipment is not in use as a regular daily occurrence – for example bain marie, display fridges etc that are turned off overnight

- Login to your dashboard
- Click on Probes Tab

This will get you to here



Probes

| Location: Castle Hill Laboratory <input type="radio"/> All <input type="radio"/> Active <input type="radio"/> Inactive | | | | | | | | | | | | | |
|--|----------------|--------------|------------|-----------|--------------|-------------|------------------|------|--------|--------|-----------|-----------|--|
| Location | Name | Serial | Alert High | Alert Low | Warning High | Warning Low | Next Calibration | Edit | Status | Action | Timeclock | Expire on | <input type="checkbox"/> Renew Online Monitoring |
| Castle Hill Laboratory | Incubator 1 | SP- BKP_0004 | 33 | 27 | 32 | 28 | 31 October 2017 | | | | | - | |
| Castle Hill Laboratory | Incubator 3 | SP- BKP_0001 | 39 | 31 | 37 | 33 | 31 October 2017 | | | | | Expired | |
| Castle Hill Laboratory | Incubator 4 | SP- BKP_0002 | 28 | 17.5 | 26 | 19 | 31 October 2017 | | | | | - | |
| Castle Hill Laboratory | Incubator 5 | SP- PIQ_0002 | 44 | 22 | 42 | 25 | 27 March 2018 | | | | | - | |
| Castle Hill Laboratory | Refrigerator 2 | SP- BKP_0010 | 9 | -2.5 | 7.5 | -1.5 | 31 October 2017 | | | | | - | |
| Castle Hill Laboratory | Refrigerator 4 | SP- BKP_0007 | 8 | -1.5 | 7 | -0.5 | 31 October 2017 | | | | | - | |

Choose data logger you want to set
Click on the 'Time Clock' icon

Choose



Warning!! Only use this feature to set data loggers that you want to turn off at certain times. Data loggers that aren't set here will run 24/7

TIME CLOCK MENU

Here is where you set the times you want your data logger to **read**.

Out of these pre-set times the data logger will be turned off.

- Choose times from drop down menu for days you want to set & press **submit**
- To stop using time clock – choose from dropdown menu START TIME 00.00 and type in END TIME 23.59

Probe Set Timeclock

[List Probes](#)

Probe Name: Incubator 1
Serial Number: SP-BKP_0004

Note: Default start time is 00:00 and end time is 23:59 which means the probe will be active for that whole day.

| Day | Start Time | End Time |
|-----------|------------------------------------|------------------------------------|
| Sunday | <input type="text" value="00:00"/> | <input type="text" value="23:59"/> |
| Monday | <input type="text" value="00:00"/> | <input type="text" value="23:59"/> |
| Tuesday | <input type="text" value="00:00"/> | <input type="text" value="23:59"/> |
| Wednesday | <input type="text" value="00:00"/> | <input type="text" value="23:59"/> |
| Thursday | <input type="text" value="00:00"/> | <input type="text" value="23:59"/> |
| Friday | <input type="text" value="00:00"/> | <input type="text" value="23:59"/> |
| Saturday | <input type="text" value="00:00"/> | <input type="text" value="23:59"/> |

← Choose



Warning!! Only use this feature to set data loggers that you want to turn off at certain times. Data loggers that aren't set here will run 24/7

ASSISTANCE

We want you to experience the full benefits of



If you need assistance or have any questions, please do not hesitate to contact us
at anytime

helpdesk@smartprobe.com.au

M 0410 557 727

O 02 8605 2942

www.smartprobe.net.au